



# Safe Place® Software Version 7.0

## Smart Sense™ Technology

A cutting edge technology that senses when the Transmitter has been removed from an infant. Adjustable alarming parameters to optimize security by unit, or by individual patient. A Red Alarm displays when the transmitter is tampered with or is loose.

## Joint Commission Reports

Two levels of 'Tracer' reports for Joint Commission. For ease-of-use, the Safe Place Reports List has been divided into sections based on the depth of reporting requested.

## Workflow/Easy-to-Use Interface

Touchscreen buttons with colorful icons for easier identification of function. Condensed menus and submenus for easy navigation.

## Pre-Enroll & Auto-Enroll Features

Pre-Enroll functionality saves time by entering patient information before infant is banded. It can be updated later with a Transmitter that is in use. Auto-Enrolls will appear on the Census Screen for accurate accounting of patients banded, and for easy identification of Safe Place Transmitters needing to be named.

## Escort Location Feature

When a patient is assigned to Escort status, the end-user selects a timeframe and destination. On the Safe Place Census Screen, the patient is reported in Escort status in one column, with the intended destination in another column.

## Multiple Alarm Notification Options

Messaging Functions including receiving Alarms via Text Message, Email, Pagers and Phone Systems including Cisco, Spectralink, Avaya, Vocera, and others.

## Event Cause Documentation

Event Resolution or Clear Reasons are configurable, allowing a facility to customize Event Reasons such as "Unauthorized Open", "Test," and "Too Close to Door".

## Login Permissions

The system recognizes different user levels, which are defined by the facility, giving staff access to specific system functions.

## System Troubleshooter

Three consecutive "No Signal" alarms launches a window with intended actions: Help, Close, Tag In-Use, Tag Not-in-Use, and Discharge; each window informs staff on how to handle the alarm.

## Equipment Monitoring/Blue Alarms

The Blue Alarm function is designed for equipment containment and location within a monitored area. Monitored equipment is also displayed on the Census Screen and are easily distinguished from patients.



Easy-to-Use Icons/Buttons



Joint Commission Reports



Pre-Enroll Feature



Transmitter Troubleshooter

RF Technologies Leads with 20 Years of Experience Designing & Delivering Innovative & Customizable RFID Solutions to the Healthcare Industry

3125 North 126th Street, Brookfield, WI 53005 USA  
phone. 800.669.9946 fax. 262.790.1784

[www.rft.com](http://www.rft.com)

# SAFE PLACE®

Infant Security Solutions