



Features

- Once activated, device remains in alarm until cleared at the point of care
- Call data can be sent to optional communication devices such as pagers, Quick Look™ LED displays and wireless phones
- Single or dual occupancy rooms
- Standard 0.25 inch plug on each device
- Large handles are easy to grasp
- Conforms to UL Standards

A variety of nurse call options are available for your CODE ALERT® Quick Response® Wireless Call Solution. Residents can quickly and easily summon help through the nurse call method that best suits their motor skills and range of movement.

- **Standard Call Cord (Momentary):** Red button clicks when pressed, assuring residents that their call has been made.
- **Standard Call Cord (Locking):** When pushed, button locks in place until reset at resident's bedside.
- **Pad Call Cord:** Activated by applying minimal pressure or leaning on a soft gray pad.
- **Touch Activated Call Cord:** Similar to a standard nurse call cord, the device is activated by the push of a rectangular pad.

Technical Specifications

Standard Nurse Call Cord - Momentary

Color/Composition	White with red button, latex free
Cord Length	10 feet
Connector	.25" phone plug
Max Voltage	28 VAC/DC
Standards	Meets UL 1069 oxygen safety requirements



Momentary Call Cord (UL)

Standard Nurse Call Cord - Locking

Color/Composition	Grey, latex free
Cord Length	10 feet
Connector	.25" phone plug
Max Voltage	28 VAC/DC
Standards	Meets UL 1069 oxygen safety requirements
Options	Clip for attaching to bed linens



Standard Momentary Call Cord



Locking Call Cord

Technical Specifications

Pad Call Cord


Color/Composition	Grey, 3" diameter, latex free
Cord Length	10 feet
Connector	.25" phone plug
Max Voltage	28 VAC/DC
Standards	Meets UL 1069 oxygen safety requirements
Options	Gripping feet on the bottom, clip for attaching to bed linens



Pad Call Cord

Touch Activated Call Cord

Color/Composition	Grey, latex free
Cord Length	10 feet
Connector	.25" phone plug
Max Voltage	28 VAC/DC
Standards	Meets UL 1069 oxygen safety requirements
Options	Clip for attaching to bed linens



Touch Activated Call Cord



Y-Adapter for dual occupancy

Selection Chart

0490-0501	Call Cord Set, Locking, Single Station, 10ft
0490-0502	Call Cord Set, Locking, Double Station, 10ft
0490-0538	Call Cord Set, Momentary, Single Station, 10ft
0490-0549	Call Cord Set, Momentary, Single Station, 10ft (UL)
0490-0505	Touch Activated Call Cord
0490-0539	Pad Call Cord, Single Station, 10 ft, Momentary
0460-0005	Nurse Call Y-Adapter

