



The SENSATEC® fall monitor is used with sensor options to alert staff if an at-risk patient or resident tries to get up without assistance. When weight is removed from a pad, the sensor is triggered and an alarm sounds at the fall monitor.

If you have a nurse call system in place, the ST750 fall monitor can connect to most bedside nurse call cord systems through a universal call jack or adapter. When a patient triggers a sensor or presses the call cord button, an alert is sent through the call system and a silent (volume set at the lowest setting), local alarm is generated.

The fall monitor can easily hang or hook on beds, chairs, or wheelchairs. The unit includes an LED on the face of the device to indicate system status.

Technical Specifications

Dimensions (LxWxD)	4.7" x 2.8" x 1.3"
Weight	0.4 lbs
Color / Composition	White monitor, blue "Hold" button
Power Requirement	(3) AA Alkaline Batteries Power Adapter, 5VDC, 500mA (optional - ST750 only)
Status Indicator	Solid Green: 30 second Hold Flashing Red: In Alarm Blinking Red (every 5 seconds): Low Battery
Jack Connections	RJ9 Modular Sensor Input Jack 0.25" Stereo Floor Sensor Jack 0.25" Stereo Nurse Call Jack (ST750 only)
Audio Output	78 dBm
Operating Temperature	41-104°F
Warranty	One year
Audible Alarm Options	Tones: 8 different Melodies: choose from 8 different (ST750 only) Voice: 10 sec voice recording (ST750 only)

Features

- "Hold" button with LED silences alarm or suspends the system and indicates status
- Begins monitoring immediately after a patient or resident is positioned on a sensor pad
- Tamper resistant: all functions are secured behind a panel with screw
- Restraint free to respect patient dignity
- Silicone protective boot
- Conforms to AAMI ES60601-1

Features	Models	
	ST710	ST750
Auto Reset	X	X
Low Battery Signal	X	X
Auto On	X	X
Adjustable Volume	X	X
Adjustable Delay	X	X
Adjustable Alarm Tone	X	X
Battery Powered	X	X
AC Power Supply		X
Voice Recording/Alert		X
Nurse Call Integration		X



(Continued)

Selection Chart

1000-0570	Sensatec ST710 Fall Monitor
1000-0575	Sensatec ST750 Fall Monitor with Voice Alert
0180-0068	AC Adaptor for ST750
0460-0150	Nurse Call Connection Cord for ST750
0900-0222	Nylon Hanging Strap, ST700 Series
0120-0162	Wall Bracket, ST700 Series
0120-0163	Mounting Clip, ST700 Series
0220-0318	Replacement Boot, ST700 Series
0460-0005	Nurse Call Y Adapter for ST750
0800-0343	Velcro Straps and Strips

