



## Features

- Helps prevent accidental infant switching
- Two-year battery life
- Receive visual match alerts before reaching the mother
- HIPAA-compliant, fulfills Joint Commission requirements

### Know Before Even Entering the Room

Today's increasingly busy obstetric departments operate with fewer staff members covering a higher workload, introducing the possibility of mismatching a newborn with his or her mother. The wireless SAFE PLACE® Baby Check® provides an added layer of protection against human error and a potential mismatch, reducing a hospital's liability and giving patients added peace of mind.

### Simple, Automatic Process Saves Time

Baby Check uses wall mounted transmitters and your existing SAFE PLACE\* transmitters to confirm the correct baby is presented to the right room and mother. Confirmation occurs before entering the room to help maintain the mother's feeling of security and confidence.

- A small transmitter is mounted outside the mother's room that corresponds with the SAFE PLACE transmitter already secured on her infant
- When a baby is brought to the door within 15 inches, nurses press the check button; a green light flashes on the Baby Check transmitter to indicate a match
- In the case of a mismatch, a red light will flash to notify the nurse

### Integrates with Existing Safe Place System

With your existing SAFE PLACE system, nurses are notified of a mismatch on a central computer station or other mobile notification devices and can respond quickly. Activities and Events are also recorded in software for automated reporting to help your hospital meet Joint Commission requirements.

*\*Compatible with Safe Place software versions 7 and higher.*

